## EDAN







## **EDAN** Edan Instruments, Inc.

- 3/F B, Nanshan Medical Equipment Park, Nanhai Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China
- Tel: +86-755-26898326 Fax; +86-755-26898330 www.edan.com.cn Email: info@edan.com.cn



## EDAN With its bright LED cold lighting system, high-resolution camera and electronic green filter, C3/C6 Video Colposcope realizes full-screen display of precise cervical examination, offering crucial and accurate information for clinical decision making to the gynecologist. ✓ Features ■ LED cold lighting system reproduces an original view of tissue Advanced and fast auto-focusing system ensures clear image Remote capture control and user-friendly keyboard layout for easy operation ■ Electronic green filter enhances vascular image display without compromise of illumination ■ Acetic acid and iodine reaction test timer to capture the right and precise image ■ Real-time display of magnification ■ Powerful image display and data management software Electronic green filter EDAN Acetic acid and iodine reaction

We also offer you patient management software especially designed for the needs of gynecologists, which enables you to take full advantage of the Video Colposcope. The software enables you to capture, review, analyze, store, print, and export colposcopic images quickly and easily.

· Patient tracking and data management

C3/C6 Video Colposcope

- · Built-in marking system for analysis
- Export report in PDF format
- Print report on A4 size paper via external printer



LED cold ligh